



MONO-VALVE CERVICAL COLLAR KIT (6 PIECES) NEC YORK KIT

FEATURES

- Nec York is **fast to apply** thanks to the embedded strap closure
- **Comfortable shaped padding**
- It provides **access for tracheotomy** and **openings in the nuchal section** for better transpiration
- Nec York kit includes **6 collars of 5 different sizes** (Infant, Paediatric, Neckless, Short x 2, Regular).



The **NEC YORK KIT** comes with 6 pieces of mono-valve cervical collars, plus transport bag. Each collar is colour-coded depending on its size (5 sizes in total). It is built to maintain a neutral alignment of the patient cervical spine.

TECHNICAL DATA¹

Kit weight with bag / One piece weight	0,76 kg / 0,12 kg
Size 1 (Infant)	12,5 x 45 x 8 cm
Size 2 (Paediatric)	12,5 x 44,5 x 10 cm
Size 3 (Neckless)	13,5 x 55,5 x 9,5 cm
Size 4 (Short)	13 x 55,5 x 8 cm
Size 5 (Regular)	14,5 x 55,5 x 11,5 cm
Materials	PE, PA

¹ Length measures are subjected to tolerances of ± 10 mm

ACCESSORIES



QC10120A

COLLAR BAG 1

Bag for mono-valve collars **[INCLUDED]**



QC10120A

COLLAR BAG 3

Bag with dividers mono-valve collars

CLASS I MEDICAL DEVICE COMPLIANT WITH UE REGULATION 2017/745

Product datasheet **NEC YORK KIT** - Code **QC19015** including: Codes **QC19001** (size 1 - Infant) • **QC19002** (size 2 - Paediatric) • **QC1902C** (size 3 - Neckless) • **QC19021** 2 x size 4 - Short) • **QC19022** (size 5 - Regular)

Rev.0 (31/03/2021)

©Copyright - This document is property of SPENCER ITALIA S.r.l. and can't be reproduced, used and disclosed without the written approval of SPENCER ITALIA SRL.

Spencer Italia s.r.l. a socio unico soggetta a direzione e coordinamento da Protect Medical Holding GmbH
Via Provinciale, n° 12 - 43038 Sala Baganza (PR) Italy - +39 0521 541111 - info@spencer.it - www.spencer.it